

# Example

December 7, 2020

## 1 Experimental setup

### SIRIUS Beamline : Experiment 123456789

**This is an example of JupyLabBook**

- Type: Proposal
- Safety: Yellow
- Date: 12/03/2020-17/03/2020
- Main proposer: Doe
- Local contact: AH
- Users (on site): John Doe, Jane Doe
- Recording directory: /Users/arnaudhemmerle/Documents/Recherche/Analysis/JupyLabBook/recording/
- Machine:
  - Current: 450 mA
  - Mode: Top-up
- Optics:
  - DCM: Si111
  - MGM: Not used
  - M1: M1-A Pt Track
  - M2: M2 Pt Track
  - M3: No M3
  - M4: M4 Pt Track
- Beam:
  - Fixed/Variable energy: Fixed
  - Energy (keV): 8
  - Wavelength (nm): 0.155
  - Harmonic: 19
  - Polarisation: LH
  - Phase (deg): 0
  - Horizontal focalisation: False
  - Vertical focalisation: True
  - Horizontal beamsize (mm): 2
  - Vertical beamsize (mm): ~0.150
- Monitors and XBPM:
  - mon1:

- mon2: thick diamond
  - mon3:
  - mon4:
  - Detectors: Pilatus on Delta-Gamma
- Remarks: This is an example.

## 2 Beamline alignment

Here we show functions used during beamline alignment.


LaTeX formula can be used:

$$\frac{786 - 558}{2 \times 2069} \times 0.0355 = 1.9 \text{ mrad}$$

### 2.1 Subsection


#### 2.1.1 SIRIUS\_2020\_03\_11\_0744: dscan basez -.2 .2 50 .1

Fit with erf function.




### 2.1.2 SIRIUS\_2020\_03\_11\_0752: One can edit the scan info here as well

Fit with a gaussian.



## 2.2 Calibration thetaz



## 3 GIXD

### 3.0.1 SIRIUS\_2020\_03\_12\_0756: continuous\_ascan delta -24 -19 100 5

Extraction of the Vineyard.

#### - Open Nexus Data File :




```
recording/SIRIUS_2020_03_12_0756.nxs
. Number of data points: 101
. Available Counters:
  0 -----> delta
  1 -----> zs
  2 -----> gamma
  3 -----> hu36energy
  4 -----> xs
  5 -----> energydcm
  6 -----> current
  7 -----> mon2
  8 -----> surfacepressure
  9 -----> areapermolecule
 10 -----> qxy
```

```

11 -----> pilatus
12 -----> pilatusroi1
13 -----> integration_time
14 -----> sensorsRelTimestamps
15 -----> sensorsTimestamps
. Pilatus data found, (column 11, alias pilatus)
. qxy data found, (column 10, alias qxy)
. Surface pressure data found, mean value 19.74 ± 0.006119 mN/m
. Area per molecule data found, mean value 0.3557 ± 3.944e-05 nm2 per
molecule
. Gamma motor data found, mean value -0.0004471 deg

```

SIRIUS\_2020\_03\_12\_0756.nxs



Data not saved. To save data, run a GIXD on the scan.  
Channel10: 607

### 3.0.2 SIRIUS\_2020\_03\_12\_0756: continuous\_ascan delta -24 -19 100 5


Classic GIXD with:

$$q_{xy} = \frac{4\pi}{\lambda} \sin\left(\frac{2\theta}{2}\right)$$

Generates:

- SIRIUS\_2020\_03\_12\_0756\_1D\_qz.dat for each binning
- SIRIUS\_2020\_03\_12\_0756\_1D.dat
- SIRIUS\_2020\_03\_12\_0756\_1D.mat for each binning
- SIRIUS\_2020\_03\_12\_0756\_1D.moy for each binning

. Absorbers: 29 - Vide



### 3.0.3 SIRIUS\_2020\_03\_12\_0756: continuous\_ascan delta -24 -19 100 5

It is possible to print all info on the scan and the counters.

- Open Nexus Data File :

recording/SIRIUS\_2020\_03\_12\_0756.nxs

. Number of data points: 101

. Available Counters:


0	----->	delta
1	----->	zs
2	----->	gamma
3	----->	hu36energy
4	----->	xs
5	----->	energydcm
6	----->	current
7	----->	mon2
8	----->	surfacepressure

```


    9 -----> areapermolecule
   10 -----> qxy
   11 -----> pilatus
   12 -----> pilatusroi1
   13 -----> integration_time
   14 -----> sensorsRelTimestamps
   15 -----> sensorsTimestamps
. Absorbers: 29 - Vide
. Pilatus data found, (column 11, alias pilatus)
. qxy data found, (column 10, alias qxy)
. Valid data between points 0 and 100
. Surface pressure data found, mean value 19.74 ± 0.006163 mN/m
. Area per molecule data found, mean value 0.3557 ± 3.866e-05 nm2 per
molecule
. Gamma motor data found, mean value -0.0004715 deg
. Original, non binned matrix saved in:
working/SIRIUS_2020_03_12_0756_1D.mat
. Scalar data saved in:
working/SIRIUS_2020_03_12_0756_1D.dat
. Qz values saved in:
working/SIRIUS_2020_03_12_0756_1D_qz.dat10
. Binned matrix saved in:
working/SIRIUS_2020_03_12_0756_1D.mat10
. XYZ    data saved in:
working/SIRIUS_2020_03_12_0756_1D.moy10
. Qz values saved in:
working/SIRIUS_2020_03_12_0756_1D_qz.dat20
. Binned matrix saved in:
working/SIRIUS_2020_03_12_0756_1D.mat20
. XYZ    data saved in:
working/SIRIUS_2020_03_12_0756_1D.moy20
. Qz values saved in:
working/SIRIUS_2020_03_12_0756_1D_qz.dat40
. Binned matrix saved in:
working/SIRIUS_2020_03_12_0756_1D.mat40
. XYZ    data saved in:
working/SIRIUS_2020_03_12_0756_1D.moy40

```

SIRIUS\_2020\_03\_12\_0756.nxs



$\pi = 19.74 \text{ mN} \cdot \text{m}^{-1}$





$\gamma = -0.0004715 \text{ deg}$

## 4 Isotherm

### 4.0.1 SIRIUS\_Isotherm\_2019\_02\_17\_01544: isotherm 1.97 46 35000 1

Plot the isotherm. Generates SIRIUS\_Isotherm\_2019\_02\_17\_01544.dat

SIRIUS\_Isotherm\_2019\_02\_17\_01544






## 5 1D plot

### 5.0.1 SIRIUS\_2020\_03\_12\_0760: run cont\_regh.ipynb

Add a 1D plot by clicking on "Add plot to report".

Generates SIRIUS\_2020\_03\_12\_0760.dat



## 6 GIXS

### 6.0.1 SIRIUS\_2019\_11\_07\_00325: tscan 10 10

GIKS:  $q_z$  vs  $q_{xy}$ .


Image and profiles with the approximation  $q_{xy} = \frac{4\pi}{\lambda} \sin\left(\frac{2\theta}{2}\right)$ .

Generates:

- SIRIUS\_2019\_11\_07\_00325\_pilatus\_sum.tiff
- SIRIUS\_2019\_11\_07\_00325\_pilatus\_sum.mat
- SIRIUS\_2019\_11\_07\_00325\_integrated\_qz.dat
- SIRIUS\_2019\_11\_07\_00325\_integrated\_qxy.dat

## 6.0.2 SIRIUS\_2019\_11\_07\_00325: tscan 10 10

```
. Absorbers: 29 - Vide  
  
. No gamma found! gamma = 0  
. No delta found! delta = 11.578
```



## 7 Plot pilatus


### 7.0.1 SIRIUS\_2019\_11\_07\_00325: tscan 10 10

Plot the sum of the images from the Pilatus (time integration).

Generates:

- SIRIUS\_2019\_11\_07\_00325\_pilatus\_sum.tiff
- SIRIUS\_2019\_11\_07\_00325\_pilatus\_sum.mat
- SIRIUS\_2019\_11\_07\_00325\_integrated\_x.dat
- SIRIUS\_2019\_11\_07\_00325\_integrated\_y.dat

```
. Absorbers: 29 - Vide
```



## 8 XRF


### 8.0.1 SIRIUS\_2017\_12\_11\_08042: run xsw7.ipynb



Plot XRF from the 4-elements detector, in channels and without peak identification.

Generates:

- SIRIUS\_2017\_12\_11\_08042\_fluospectrum.mat for each element
- SIRIUS\_2017\_12\_11\_08042.dat

. Absorbers: Al 200micron








### 8.0.2 SIRIUS\_Fluo\_2020\_07\_03\_0042: tscan 500 30

Plot XRF from the 1-element detector, in eVs and with peak identification.

Generates:

- SIRIUS\_Fluo\_2020\_07\_03\_0042\_fluospectrum04.mat
- SIRIUS\_Fluo\_2020\_07\_03\_0042.dat

. Absorbers: Al 800micron



## 9 Insert script

Script inserted (with automatic scan numbering) using “Insert script”.

### 9.0.1 script\_with\_loop.ipynb

```
%shopen
%amove delta -40
%run reset_motors.ipynb
%amove delta -35
%continuous_ascan delta -35 -25 250 5 #123
%run reset_motors.ipynb
for i in range(4):
    %amove delta -20
    %continuous_ascan delta -10 -3 175 5 #124 #126 #128 #130
    %run reset_motors.ipynb
    %run cont_regh_abs.ipynb #125 #127 #129 #131
for i in range(3):
    %amove delta -20
    %continuous_ascan delta -10 -3 175 5 #132 #135 #138
    %run reset_motors.ipynb
    %run cont_regh_abs.ipynb #133 #136 #139
    %run cont_regh_abs.ipynb #134 #137 #140
for i in range(2):
    %amove delta -20
    %tscan 10 100 #141 #142

%continuous_ascan delta -35 -25 250 5 #143
%amove delta -40
%shclose
```

## 10 Insert positions

Positions extracted from the logs, using “Insert positions”.

### 10.0.1 wm zs

zs
-40.9178

### 10.0.2 wm diffracto

deltacodeur	euchi	euth	euphi	kappa_h	kappa_k
	1.00196	-89.57961	90.42039	-0.00580	-0.08252
Degrees	deg	deg	deg		

kappa_l	qxy	qxy0	qz	basexPoint	basexTrait
-0.18486	2.0556	23.82	-0.92	-15.7275	-15.7274
	nm-1	nm-1	nm-1		

basezPlan	basezPoint	basezTrait	basepitch	baseroll	basex
71.1257	71.1257	71.1257	-0.0000	0.000	-15.727
			mmrad	mmrad	mm

baseyaw	basez	alphax	alphay	delta	delta0
-0.000	71.126	0.2998	0.2000	-2.9110	-34.2322
mmrad	mm				

deltaa	etaa	gamma	kappav	mu	kphi
0.0000	0.0000	1.2997	1.3080	-179.9997	0.0000

thetaa	thetah	komega	xs	ky	ys
0.0000	0.0185	0.0000	0.0000	-0.1000	0.0000

kz	zs	kx
0.0000	-41.9999	-0.1000

## 11 Insert commands

Commands extracted from the logs, using "Insert commands".

```

Wed, 11 Mar 2020 13:59:25 dmove m4tz 1
Wed, 11 Mar 2020 13:59:29 ct 1
Wed, 11 Mar 2020 13:59:37 amove m4pitch -.02
Wed, 11 Mar 2020 13:59:39 ct 1

```

```

Wed, 11 Mar 2020 13:59:50 dmove m4tz -.5
Wed, 11 Mar 2020 13:59:53 ct 1
Wed, 11 Mar 2020 13:59:56 dmove m4tz -.5
Wed, 11 Mar 2020 14:00:00 ct 1
Wed, 11 Mar 2020 14:00:12 dscan m4tz -.5 1.5 50 .1 #679
Wed, 11 Mar 2020 14:01:47 sigmoid_dscan m4tz -.4 .4 50 .1 #680
Wed, 11 Mar 2020 14:03:41 amove m4tz PEAK
Wed, 11 Mar 2020 14:03:50 dscan m4pitch -.5 .5 50 .1 #681
Wed, 11 Mar 2020 14:05:04 sigmoid_dscan m4tz -.4 .4 50 .1 #682


```

## 12 Convert logs

Human-readable logs generated in the folder `/working/readable_logs/` by clicking on “Convert logs”.

## 13 Insert an image

Using the command “Insert image”.





## 14 Export to pdf

PDF generated by clicking on ''Export to pdf''.